

# ALTAPURE<sup>®</sup>



## Making Things Clean: Saving Lives

### Introducing the **HLDC-1**

### A Multi-Purpose High Level Disinfection Cabinet

The HLDC-1™ is an advanced multi - purpose high-level disinfection cabinet, using Altapure's enhanced ultrasonic sub-micron aerosol technology, delivering a consistent 100% kill of pathogens such as: C. difficile spores, VRE, CRE, MRSA, C. auris, & viruses such as HPV.

#### FEATURES:

- Rapid Process: less than a 20 minute cycle time.
- Economical: less than \$0.20 per cycle.
- Efficacious: delivers a 100% kill of C. diff spores, MRSA, VRE, CRE, HPV, and C. auris.
- Gentle On Equipment: very low concentration of PAA, low temperature / cool process.
- Easy To Use: computer controlled.
- Data: Wireless data reporting.
- Specifications: average droplet size 0.69 micron.
- Electrical: 120VAC 50/60Hz 15A. All global standard input options are available.
- Dimensions: 1. Outside: 72 X 36 X 24 inch  
2. Chamber: 54 X 27.7 X 12 inch
- Construction: all food grade materials.
- Disinfects: entire ultrasound probes, patient monitoring cables, stethoscopes, key boards, phones, blood pressure cuffs, and much more.



## IMPROVING THE PATIENT EXPERIENCE